



Total Blood Bank Processing Solutions

PLASMA SOFT

DESCONGELADOR Y CALENTADOR DE SANGRE - PLASMA THAWING

El descongelador y calentador de sangre PLASMA SOFT permite el descongelamiento de plasma en forma controlada y segura, y el calentamiento de hemoderivados para llevarlos a temperatura de transfusión.

PLASMA SOFT is a revolutionary thawing system, which provides controlled and safe plasma defrosting, together with the heating of any blood products that need to be taken to an optimum temperature in order to be transfused.

[1e] 28 38 48
58 68 78 88

BOLSA Min:Seg
2 '12': 10

[1e] 28 38 48
58 68 78 88

12 [28] 38 48
58 68 78 88

SET. TEMPERATURA
'37.0 °C'



PLASMA SOFT

Un software especialmente diseñado, combinado con un sistema de transferencia de calor hacia la unidad, permite el descongelamiento y calentamiento individual, brindando seguridad y control en el proceso.

To achieve maximum safety and control of the process, we have developed an exclusive software that works together with an individual heating system, enabling the thawing and heating of each unit to be done separately.



Calle 34 (Ex Francia) N° 3917 - San Martín - B1650IHM - Buenos Aires - Argentina

Tel/Fax: (54-11) 4755-7780/ 4755-4004/ 4754-3532

http://www.presvac.net - e-mail: presvac@presvac.net

