

SILLON ARTICULADO ELÉCTRICO - ELECTRIC ADJUSTABLE CHAIR



T-5 E

La línea de sillones de extracción Presvac, tiene todas las características técnicas necesarias para la extracción de sangre del donante, de manera confortable y segura, además de brindarle al operador una óptima y eficiente postura de trabajo.

- Tapicería de cuero ecológico de primera calidad
- Ajuste eléctrico de los ángulos
- Posición Trendelenburg
- Funcionalidad de RCP
- Pedal, reposacabezas y portasueros
- Estructura de acero con recubrimiento en polvo
- Ruedas de 4 pulgadas
- Diseño lujoso y cómodo

Material	Acero pintado PVC
Dimensiones (abierto)	190×88×60cm
Ruedas	4" ruedas con freno
Accesorios	soporte para suero, reposacabezas y escalón para pies
Límite de peso	200KG
Dimensiones embalaje	1pcs/Carton, 151×82×60cm
N.W.	62kg
G.W.	72kg

The Presvac blood donation chair line offers all the technical features required for safe and comfortable blood collection, while also providing the operator with an optimal and efficient working posture.

- Premium eco-friendly leather upholstery
- Durable powder-coated steel frame
- Electrically adjustable angles
- Trendelenburg function
- 4-inch casters with brakes for easy mobility
- Equipped with footrest, headrest, and IV pole
- Luxurious design and superior comfort



Relaxing Position



Semi-Fowler Position



Lying Position



Trendelenburg Position (~12°)

